



The Bethesda System 2001 Reporting System for Cervical Cytology

Stuart Dobson
Manager, Customer Education & Support



HOLOGIC CLARITY OF VISION

HOLOGIC CLARITY OF VISION

Why use The Bethesda System 2001?

- A single global reporting system facilitates international epidemiological research
- We now have greater understanding of the pre-invasive phase of cervical cancer.
 - LSIL(CIN1/HPV) = Acute viral infection
 - HSIL = HPV induced cx. cancer precursor
- TBS 2001 diagnostic categories better reflect patient management options
 - ASC-US = Repeat cytology or HPV testing
 - LSIL = Repeat cytology 6 months
 - HSIL = Colposcopy and Biopsy
- Reduced diagnostic categories improves inter- and intra-observer variability



HOLOGIC CLARITY OF VISION

Comparison of Reporting Nomenclature for Cervical Cytology

Papanicolaou System	Class I, II	Class IIIa		Class IIIb	Class IV	Class V	
	ASC-US	ASC-H					
WHO	Negative	Mild Dysplasia (CIN 1)	Moderate Dysplasia (CIN 2)	Severe Dysplasia (CIN 3)	CIS	Microinv. SCC	Invasive SCC
TBS 2001	Negative	LSIL	HSIL			HSIL (susp. for invasion)	SCC
		ASC-US	ASC-H				
		AGC(NOS). AGC(favor neoplastic)			AIS	Adenocarcinoma	



The Bethesda System 2001

SPECIMEN TYPE

Indicate conventional smear (Pap smear) vs. liquid based vs. other


SPECIMEN ADEQUACY

- Satisfactory for evaluation
(describe presence or absence of transformation zone component and any other quality indicators, e.g., partially obscuring blood, inflammation, etc.)
- Unsatisfactory for evaluation ... (specify reason)



Conventional Pap Smear Specimen Adequacy

- Approx. 8,000 – 12,000 well preserved squamous cells
- An estimation based on reference images below



ThinPrep® Pap Test Specimen Adequacy

The Bethesda System 2001 Cellular Composition for Adequacy

TP Diam mm	Area mm ²	FN20 eyepiece/ 10x obj.		FN20 eyepiece/ 40x obj.		FN22 eyepiece/ 10x obj.		FN22 eyepiece/ 40x obj.	
		# Fields @ FN20 10X	# Cells/Field for 5K Total	# Fields @ FN20 40X	# Cells/Field for 5K Total	# Fields @ FN22 10X	# Cells/Field for 5K Total	# Fields @ FN22 40X	# Cells/Field for 5K Total
20	314.2	100	50.0	1600	3.1	82.6	60.5	1322	3.8

- A minimum of 5000 well preserved squamous cells



The Bethesda System 2001

NEGATIVE FOR INTRAEPITHELIAL LESION OR MALIGNANCY (NILM)

- ORGANISMS:
 - *Trichomonas vaginalis*
 - Fungal organisms morphologically consistent with *Candida* spp
 - Shift in flora suggestive of bacterial vaginosis
 - Bacteria morphologically consistent with *Actinomyces* spp.
 - Cellular changes consistent with Herpes simplex virus



The Bethesda System 2001

NEGATIVE FOR INTRAEPITHELIAL LESION OR MALIGNANCY (NILM)

- OTHER NON NEOPLASTIC FINDINGS
 - Reactive cellular changes associated with
 - inflammation (includes typical repair)
 - radiation
 - intrauterine contraceptive device (IUD)
 - Glandular cells status post hysterectomy
 - Atrophy
- OTHER
 - Endometrial cells (> 40 years of age)



The Bethesda System 2001

EPITHELIAL CELL ABNORMALITIES

SQUAMOUS CELL

- Atypical squamous cells
 - of undetermined significance (ASC-US)
 - cannot exclude HSIL (ASC-H)
- Low grade squamous intraepithelial lesion (LSIL) encompassing: HPV/CIN 1
- High grade squamous intraepithelial lesion (HSIL) encompassing: CIN 2 and CIN 3/CIS
 - with features suspicious for invasion
- Squamous cell carcinoma



The Bethesda System 2001

Atypical Squamous Cells of Undetermined Significance (ASC-US)

- Atypical cells falling short of LSIL (CIN 1, HPV)
- USA median laboratory reporting rate 4.3% ¹
- Conventional Pap Range 0.3% - 8.2%
- Management
 - Reflex HPV Testing if using LBC
 - Repeat Cervical Cytology in 6 months
- Risk of HSIL on biopsy = 10% ²

1. Eversole GM et al. Arch Path Lab Med, Mar 2010
2. Wright T et al, NCCN 2002



The Bethesda System 2001

Low-Grade Squamous Intraepithelial Lesion (LSIL)

- Incorporates diagnoses HPV and CIN1, mild dysplasia
- USA median laboratory reporting rate 2.5% ¹
- Conventional Pap Range 0.3% - 6.7%
- Management
 - Repeat Cervical Cytology in 6, 12 months or colposcopy
- Risk of HSIL on biopsy = 15%-30% ²

1. Eversole GM et al. Arch Path Lab Med, Mar 2010
2. Moscicki AB et al. JAMA, 2001



The Bethesda System 2001

Atypical Squamous Cells - cannot exclude HSIL (ASC-H)

- Atypical cells present where HSIL cannot be excluded
 - Cells usually scant
- USA median laboratory reporting rate 0.3% ¹
- Conventional Pap Range 0.1% - 2.0%
- Management
 - Colposcopy and Biopsy
- Risk of HSIL on Biopsy = 24 - 94% ²

1. Davey D et al. Arch Path Lab Med, Nov 2004
2. Wright T et al, NCCN 2002



The Bethesda System 2001

High-Grade Squamous Intraepithelial Lesion (HSIL)

- Incorporates diagnoses CIN2 (moderate dysplasia) and CIN3 (severe dysplasia)
- USA median laboratory reporting rate 0.5% ¹
- Conventional Pap Range 0.1% - 2.0%
- Management
 - Colposcopy and Biopsy
- Risk of HSIL on biopsy = 92-98% ²

1. Eversole GM et al. Arch Path Lab Med, Mar 2010
2. Belinson J et al. Gyn Onc 2001



The Bethesda System 2001

Squamous Cell Carcinoma

- The conventional Pap smear has proven very effective at reducing incidence of Squamous Cell Carcinoma
- Many pathologists utilise verbiage such as "*High grade epithelial abnormality, cannot exclude Squamous Cell Carcinoma*"



The Bethesda System 2001

EPITHELIAL CELL ABNORMALITIES

GLANDULAR CELL

- Atypical
 - endocervical cells (exceeds reactive, reparative)
 - endometrial cells
 - glandular cells
- Atypical
 - endocervical cells, favor neoplastic (just short of AIS)
 - glandular cells, favor neoplastic
- Endocervical adenocarcinoma in situ (AIS)
- Adenocarcinoma
 - endocervical
 - endometrial
 - extrauterine



The Bethesda System 2001

Atypical Glandular Cells (AGC)

- Endocervical
 - Endometrial
 - Glandular
- } Exceeds diagnosis of reactive, reparative change
- USA median reporting rate 0.2% ¹
 - Management
 - Colposcopy plus endocervical sampling
 - Women >35 (with abn. bleeding or atypical endom.) = colposcopy plus endocervical & endometrial sampling
 - Risk of HSIL, AIS, or Ca on Biopsy = 9 - 41% ²

1. Davey D et al. Arch Path Lab Med, Nov 2004
2. Wright T et al, NCCN 2002



The Bethesda System 2001

Atypical Glandular Cells – Endocervical

- May reflect the presence of reparative/regenerative process or cervicitis
- Typically used for cases of repair with atypia

Atypical Glandular Cells – Endocervical

- May reflect the presence of endometrial hyperplasia or proliferative endometrium
- Typically used when endometrial cells are present in peri- or post-menopausal patients not on HRT
- Cells may appear cytologically benign



The Bethesda System 2001

Atypical Glandular Cells – favour neoplastic

- **Endocervical** } Falling short of definitive diagnosis of AIS or neoplasia
- **Glandular** }
- Total AGC USA median laboratory reporting rate 0.2%¹
- Management
 - Colposcopy and Endocervical Sampling, Cone Biopsy (cold knife)
- Risk of HSIL, AIS, or Ca on Biopsy = 27-96%²

1. Davey D et al. Arch Path Lab Med, Nov 2004
2. Wright T et al, NCCN 2002



The Bethesda System 2001

Adenocarcinoma in situ (AIS)

- AIS median laboratory reporting rate 0.01%¹
- Not easily identified²
- Management
 - Colposcopy and Endocervical Sampling, Cone Biopsy (cold knife)

1. Davey D et al. Arch Path Lab Med, Nov 2004
2. Renshaw A. Arch Path Lab Med, Feb 2004



The Bethesda System 2001

Endocervical Adenocarcinoma

- Now account for 15-20% of invasive cervical cancer
- Endocervical adenocarcinoma have increased 4 fold since the 1970's

Endometrial Adenocarcinoma

- Incidental finding on Pap test
- Improved sensitivity and specificity with the ThinPrep Pap Test



Summary

- **Recommendations for Specimen Adequacy**
- **Precancerous lesion division reflect current understanding of HPV infection**
LSIL, HSIL
- **Recommendation for Atypical cells**
ASC-US, ASC-H, AGC