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

Centrifuge

PrepStain™

Vortex, Computer, Tube Vacuum, Vacuum Pump

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INTENDED USE

The ThinPrep® 2000 System is intended as a replacement for the conventional method of Pap smear preparation for use in screening for the presence of atypical cells, cervical cancer, or its precursor lesions (Low-grade Squamous Intraepithelial Lesions, High-grade Squamous Intraepithelial Lesions, and Cervical Intraepithelial Neoplasia) as defined by The Bethesda System for Reporting Cervical/Vaginal Cytologic Diagnoses.^{1, 2, 3, 4}

INTENDED USE

The PrepStain™ System (formerly the AutoCyte® PREP System) is a liquid-based thin layer cell preparation process. The PrepStain™ System produces SurePath™ slides that are intended as replacements for conventional gynecologic Pap smears. SurePath™ slides (formerly AutoCyte® PREP slides) are intended for use in the screening and detection of cervical cancer, pre-cancerous lesions, atypical cells and all other cytologic categories as defined by The Bethesda System for Reporting Cervical/Vaginal Cytologic Diagnoses.^{1, 2, 3, 4}

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CONCLUSIONS

The ThinPrep® 2000 System is as effective as the conventional Pap smear in a variety of patient populations and may be used as a replacement for the conventional Pap smear method for the detection of atypical cells, cervical cancer, or its precursor lesions, as well as all other cytologic categories as defined by The Bethesda System.

The ThinPrep 2000 System is significantly more effective than the conventional Pap smear for the detection of Low-grade Squamous Intraepithelial Lesions and severe lesions in a variety of patient populations.

Specimen quality with the ThinPrep 2000 System is significantly improved over that of conventional Pap smear preparation in a variety of patient populations.


SurePath™ slides produced by the PrepStain™ System provide similar results to conventional Pap smears in split-sample comparisons in a variety of patient populations and laboratory settings. In addition, there were significantly fewer Unsatisfactory and SBCLB cases with SurePath™ slides as compared to conventional Pap smears. The SurePath™ slide may thus be used as a replacement for the conventional Pap smear for the detection of atypical cells, pre-cancerous lesions, cervical cancer, and all other cytologic categories defined by The Bethesda System.

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38

Comparison of cytological and histological (biopsy) diagnosis concordant rate to Liqui-PREP with Pap-Smear

Histological diagnosis (biopsy)	Cytological diagnosis Liqui-PREP						Total
	ASC-H	LSIL	HSIL	AGC	GCC	SCC	
Total	98	567	183	13	6	57	924
Concordant	92	541	165	8	5	54	865
% of agreement	93.9	95.4	90.2	61.5	83.3	94.7	93.6

Histological diagnosis (biopsy)	Cytological diagnosis Pap-Smear						Total
	ASC-H	LSIL	HSIL	AGC	GCC	SCC	
Total	19	128	88	5	5	19	264
Concordant	15	101	67	3	4	17	207
% of agreement	80	78.9	76.1	60	80	89.5	78.4

ASC-H-High-grade squamous intraepithelial lesion; LSIL-Low-grade squamous intraepithelial lesions; HSIL - High-grade squamous intraepithelial lesions; AGC-Atypical glandular cells; GCC-Glandular carcinoma; SCC-Squamous cell carcinoma

Nguồn: [Hao Deshou](#), Clinical utility of Liqui-PREP™ cytology system for primary cervical cancer screening in a large urban hospital setting in China, 2009 .

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Nguồn: Gayer James et al. Evaluation of Liqui-PREP™ encapsulation method for liquid-based cytology: cell loss estimates during processing. 2004.

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Table 1: Comparison of morphologic parameters

	Cellularity	Clean background	Uniform distribution	Artifacts	Cellular overlapping	Architectural	Cellular morphologic change
Liqui-PREP	Adequate	Clean	No	Rare	Yes	No	No
Pap-smear	2/3 patients required multiple samples	Not clean	Highly presented	Highly presented	Yes	Highly presented	More commonly

	%ASC-US	%ASC-H	%LSIL	%HSIL	AGC
LIQUI-PREP	6,8	0,31	1,8	0,58	0,41
Conventional Smear	2,3	0,06	0,41	0,28	0,06

Nguồn: Hao Deshou, Clinical utility of Liqui-PREP™ cytology system for primary cervical cancer screening in a large urban hospital setting in China, 2009.

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